

# **NOSE POKE RESPONSE WITH THREE COLOR CUE**

ENV-114M

USER'S MANUAL

**DOC-053**

**Rev. 1.2**

Copyright ©2012  
All Rights Reserved

Med Associates Inc.  
P.O. Box 319  
St. Albans, Vermont 05478

Phone: 802.527.2343  
Fax: 802.527.5095  
[www.med-associates.com](http://www.med-associates.com)

notes

**Table of Contents**

**Chapter 1 | Introduction ..... 1**

**Chapter 2 | CONNECTIONS..... 1**

Light Conductor Cable..... 1

Nose Poke Response Cable ..... 1

**Appendix A | Contact Information ..... 1**

## CHAPTER 1 | INTRODUCTION

This alternative response device has a 1” diameter hole that is bisected by a sensitive infrared detector placed ¼” in from the front edge. The control amplifier operates on 28 VDC and produces a ground output whenever the beam is broken. The back of the unit contains a stimulus display made up of three colored LEDs.

## CHAPTER 2 | CONNECTIONS


### Light Conductor Cable

Molex Pin #	Connector 1		Connector 2		Connector 3	
	Function	Color	Function	Color	Function	Color
1	Not Used		Not Used		Not Used	
2	Red Light Activate	Black	Yellow Light Activate	White	Green Light Activate	Green
3	+ 28 Volts	Red	Not Used		Not Used	

These connectors go to any available Output on a Standard MED Connection Panel.

### Nose Poke Response Cable

Molex Pin #	Function	Wire Color
1	28 Volt Ground	Black
2	Response (to Input Bit)	White
3	+ 28 Volts	Red



This connector goes to any available Input on a Standard MED Connection Panel.

---

**NOTE:** Refer to ENV-254 Instructions for more information on the Nose Poke Response.

---

## Appendix A | Contact Information

Please contact MED Associates, Inc. for information regarding any of our products.

Visit our website at [www.med-associates.com](http://www.med-associates.com) for contact information.

For technical questions, email [support@med-associates.com](mailto:support@med-associates.com).