

# **NOSE POKE RESPONSE**

ENV-114AM

WIRING GUIDE



**DOC-177**

**Rev. 1.1**

Copyright ©2012  
All Rights Reserved

Med Associates Inc.  
P.O. Box 319  
St. Albans, Vermont 05478

Phone: 802.527.2343  
Fax: 802.527.5095  
[www.med-associates.com](http://www.med-associates.com)


## CHAPTER 1 | OPERATION

This response device has a 1" diameter hole that is bisected by a sensitive infrared detector placed ¼" in from the front edge. The control amplifier operates on 28 VDC and produces a ground output whenever the beam is broken.

### Nose Poke Response Cable

This cable connects to any available Input on a Standard MED Connection Panel.

Molex Pin #	Function	Wire Color
1	28V Ground	Black
2	Response	White
3	+28 Volts	Red



### Dimensions

- Panel: 3.0" W x 4.875" H (7.6 cm x 12.4 cm)
- Access Opening: 1.0" x 0.4" Deep (2.5 cm x 1.0 cm)
- IR Detector Placement: 0.25" (0.64 cm) from front edge of access opening.

## APPENDIX A | CONTACT INFORMATION

Please contact MED Associates, Inc. for information regarding any of our products. Visit our website at [www.med-associates.com](http://www.med-associates.com) for contact information. For technical questions, email [support@med-associates.com](mailto:support@med-associates.com).