

MOUSE RETRACTABLE SIPPER

ENV-352AW

USER'S MANUAL

DOC-208

Rev. 1.3

Copyright ©2014
All Rights Reserved

Med Associates Inc.
P.O. Box 319
St. Albans, Vermont 05478


Phone: 802.527.2343
Fax: 802.527.5095
www.med-associates.com

HARDWARE

Control Cable

Connect the ENV-352AW 3-pin Molex Control cable to any available Output on a Standard MED Connection Panel. When the output is turned on the sipper tube extends into the chamber and remains there for as long as that output is on. Once the output is turned off, the sipper tube will retract out of the chamber.

Molex Pin #	Function	Wire Color
1	28V Ground	Black
2	MED Operate	White
3	+28 Volts	Red



Lickometer Cable

If the sipper tube is being used with an optional MED Contact Lickometer Controller, then the red lickometer cable should be connected to the MED Lickometer Controller. Refer to the User's Guide for the specific Lickometer Controller being used for further information regarding the connection of the Lickometer cable. If the sipper tube is not being used with this optional equipment, this red cable is not used.

Priming the Tube

It is important to prime the Lixit valve prior to operating. To do this, press on the tip of the valve until liquid drips out.

CONTACT INFORMATION

Please contact Med Associates, Inc. for information regarding any of our products.

Visit our website at www.med-associates.com for contact information.

For technical questions, email support@med-associates.com.