


# ILLUMINATED NOSE POKE RESPONSE

WITH OLFACTORY PORTS

ENV-375-NPP

USER'S MANUAL



**DOC-272**

**Rev. 1.0**

Copyright ©2013  
All Rights Reserved

Med Associates Inc.  
P.O. Box 319  
St. Albans, Vermont 05478

Phone: 802.527.2343  
Fax: 802.527.5095  
[www.med-associates.com](http://www.med-associates.com)

notes

**Table of Contents**

**Chapter 1 | Introduction.....1**  
Specifications ..... 1  
**Chapter 2 | Wiring.....2**  
**Chapter 3 | Operation.....3**  
Photo Beam Operation ..... 3  
**Chapter 4 | Cleaning .....4**  
**Appendix A | Contact Information.....4**

## CHAPTER 1 | INTRODUCTION

The ENV-375W-NPP and ENV-375A-NPP Illuminated Nose Poke Response with Olfactory Ports for Mouse uses an Infrared (IR) transmitter and receivers (photo beams) to detect the presence of an animal when it pokes its nose into the receptacle.

The receptacle is illuminated with a yellow LED. An output signal is sent to the device to turn the light on or off. The light is designed to be used as a cue for when a scent is detectable in the receptacle.

Olfactory ports are used to deliver and evacuate a scent stimulus.

### Specifications

#### Dimensions

ENV-375A-NPP (fits the ENV-307A chamber)

Panel: 1.7" x 3.5" (4.318cm x 8.89cm)

ENV-375W-NPP (fits the ENV-307W chamber)

Panel: 2.4" x 3.5" (6.096cm x 8.89 cm)

Cable Length: 28"

Port diameter: 3/16" barbed (tubing not included)

*Figure 1-1 ENV-375W-NPP for the wide chamber (ENV-307W)*

---




## CHAPTER 2 | WIRING

### IR Detection

Connect the 3-pin Molex connector with black, red, and white wires to any available **Input** port on a standard Med Associates Connection Panel (Smart Control or Passive). An input is generated when the IR beam is broken. The response connectors provide voltage for the light stimulus and therefore must be connected for the light to function.

*Figure 2-1 – Nose Poke Input Connection*

Molex Pin #	Function	Wire Color
1	28 Volt Ground	Black
2	Nose Poke Response	White
3	+28 Volts	Red



### Operating Specifications


Operating Voltage: 28 Volts DC  
 Operating Current: 15 ma – quiescent  
 Output (white lead): 0 Volts @ 50mA – beam broken  
 27 Volts @ 40 mA – beam unobstructed

### Light Conductor Connector

Connect the 3-Pin Molex connector with green and orange wires to any available **Output** port on the Connection Panel (eg.SG-716B). When the Output is activated, the cue light will turn on.

*Figure 2-2 – Nose Poke Output Connection*

Molex Pin #	Function	Wire Color
1	N/A	N/A
2	Operate	Green
3	+28 Volts	Orange




### CHAPTER 3 | OPERATION

The Illuminated Nose Poke Response is designed for olfactory stimulation. The olfactory stimulus is supplied via a port in the bottom of the nose poke cavity and is removed from the nose poke cavity via an exhaust port located in the roof of the cavity.

#### Photo Beam Operation

The ENV-375-NPP Illuminated Nose Poke Response utilizes an infrared photo beam for monitoring when the animal's nose enters the receptacle. When this happens, the beam is broken and a signal is sent to an input. It is important to keep the photo beam holes clear of debris so that there is no beam interference. Bright lights can affect the performance of the photo beams.

*Figure 3-1 - ENV-375-NPP full assembly and cut-away view of photo beam sensors*



## **CHAPTER 4 | CLEANING**

The nose poke receptacle should be wiped down periodically to keep it clear of debris. The photo beam holes should be cleaned with forced air.

## **APPENDIX A | CONTACT INFORMATION**

Please contact MED Associates, Inc. for information regarding any of our products.

Visit our website at [www.med-associates.com](http://www.med-associates.com) for contact information.

For technical questions, email [support@med-associates.com](mailto:support@med-associates.com).