




## SMALL ANIMAL CYSTOMETRY

Systems Overview . . . . .	2-3
System Packages . . . . .	4-5
System Components . . . . .	6-9
Interface Packages . . . . .	10-11
Interface Components . . . . .	12-13
Software . . . . .	14-15





The small animal cystometry lab station is designed for performing in-vivo measurements of bladder function in mice and rats. It is ideal for physiological assessment of bladder function in “knock-out” and transgenic mice. This system allows standard urological tests and cystometrograms to be performed in conscious, unrestrained mice and rats while measuring bladder capacity, filling pressure, micturition pressure, and voided urine volume. The experimental animal is chronically instrumented with an intra-bladder catheter, which is routed underneath the skin and exteriorized at the base of the animal’s neck.\*


A computer-controlled syringe pump infuses saline into the bladder, and an in-line pressure transducer records intra-bladder pressure during the infusion. When bladder capacity is reached, the transducer records pressure during the urination event (micturition). An analytical balance is located beneath the animal cage, which records the weight of urine to provide an index of voided urine volume.

\* This technique requires familiarity with small animal survival surgery.\*

PACKAGE CONTENTS				
A = CAT-CYT-M			B = CAT-CYT-R	
#	SKU	DESCRIPTION	A	B
1	ENV-016M-CYT	MDF Sound Attenuating Cubicle	A	B
1	CAT-ENV-272-KIT	Cystometry Cage w/Connectors for Rat		B
1	CAT-ENV-273-KIT	Cystometry Cage w/Connectors for Mouse	A	
1	ENV-026	Switch Box for Cubicle Fan and Light(s)	A	B
1	SCALE-TOR-AD-520-USB-URO	AD520 Scale w/USB Kit	A	B
1	SYR-PUMP-USB	USB Infusion Pump Kit	A	B
1	TDR-100-1	Pressure Transducer (reusable/disposable)	A	B
1	CANL-430A	Bridge Amplifier Card	A	B
1	IO-206	Direct Input Module	A	B
1	CAB-TDR-100-1	Transducer Extension Cable	A	B
1	SG-210CP-25	Power Cable, 25' (7.6 m)	A	B
1	SG-216A	3-Pin Mini-Molex Extension, 18" (45.7 cm)	A	B



PACKAGE CONTENTS				
A = CYT-M-FIT-KIT			B = CYT-R-FIT-KIT	
#	SKU	DESCRIPTION	A	B
10	FAB-CYT-102	Catheter for Rat		B
10	FAB-CYT-107	Catheter for Mouse	A	
10	HAR-CYT-106	Luer 3-Way Stopcock	A	B
10	HAR-CYT-109	Male Luer Slip Coupler	A	B
10	HAR-LUER-STUB-18G	Needle Stub, 18 Gauge	A	A
10	HAR-LUER-STUB-23G	Needle Stub, 23 Gauge	A	B
10	HAR-LUER-STUB-30G	Needle Stub, 30 Gauge	A	
5	HAR-SYRINGE-10CC	Syringe, 10 mL (cc)	A	
5	HAR-SYRINGE-20CC	Syringe, 20 mL (cc)		B
2	STO-TUBING-TYGON-.040	Tygon Tubing (0.04"ID)	A	B



CYSTOMETRY LAB STATION

CAT-CYT-M **MOUSE**  
CAT-CYT-R **RAT**

SUMMARY OF CONTENTS


- > Animal cage w/modified top
- > Kit with catheters & connectors
- > Analytical balance
- > Razel™ infusion pump
- > Pressure transducer
- > Specialized sound attenuating cubicle
- > Electronic interface with push button event markers, pressure transducer amplifier, and manometer for calibration

TUBING + FITTING KIT

CYT-M-FIT-KIT **MOUSE**  
CYT-R-FIT-KIT **RAT**

- Included in the cystometry cage kits
- Precision-made catheters with a flared tip for easy implantation into a bladder






**ENV-016M-CYT**

**WORKING AREA (w×h×d)**  
22" x 27.5" x 14"  
(55.88 x 69.85 x 35.56 cm)

**SOUND ATTENUATING CUBICLE**

**ENV-016M-CYT** **MOUSE+RAT**

- Built-in structures make for easy placement of syringe pump, cage, and scale
- Provides isolation of your animal from environment
- Windows for convenient viewing



**WORKING AREA (w×h×d)**

**ENV-272** 5.8" x 4.2" x 8.8"  
(12.7 x 10.67 x 22.35 cm)

**ENV-273** 3.3" x 2.52" x 5.2"  
(8.38 x 6.4 x 13.21 cm)

**CYSTOMETRY CAGE**

**ENV-272** **RAT** | **MESH FLOOR**  
**ENV-273** **MOUSE** | **MESH FLOOR**  
**ENV-273-GR** **MOUSE** | **GRID FLOOR**

**CYSTOMETRY CAGE w/FOOD+WATER**

**ENV-272A** **RAT** | **MESH FLOOR**  
**ENV-273A** **MOUSE** | **MESH FLOOR**  
**ENV-273A-GR** **MOUSE** | **GRID FLOOR**

- Modified top for routing catheter to the animal
- Cages for both rat and mouse fit in the same SAC without modification
- Mouse models are also available with a grid floor
- Food and water models are modified to accommodate an ad-lib food cup and a volumetric drinking tube with Lixit valve



**ANALOG OUTPUT MODULE**  
IO-200-4 **MOUSE+RAT**

Four BNC ports take analog signals from your cystometry system and output them to a third-party data collection system.

Especially useful for recording bladder pressure simultaneously with other physiological signals, such as EMG or EEG electrophysiological recordings.

- Channels available for output:
  - Bladder pressure
  - Event marks

**PRESSURE TRANSDUCER SIGNAL CONDITIONER MODULE**  
CANL-430A **MOUSE+RAT**

- Used for recording bladder pressure
- A second module can be added for recording of an additional pressure (blood, abdominal, etc.)

**PUSH BUTTON DIRECT INPUT MODULE**  
IO-206 **MOUSE+RAT**

NOTE: All of the following functions are also able to be executed via MED-CMG software.

Four push buttons to manually:

- Start the infusion pump
- Stop the infusion pump
- Annotate a micturition event
  - Infused volume trace is reset to zero and counting begins again for the next micturition cycle
- Annotate a recording artifact
  - Adds an event mark to the data that alerts the user that something other than urine has disturbed the scale



SCALE-TOR-AD-520-USB-URO

SYR-PUMP-USB

**ANALYTICAL BALANCE**  
SCALE-TOR-AD-520-USB-URO **MOUSE+RAT**

Analytical balance for recording voided urine volume. Contains a USB cable for easy connection to data acquisition computer.

**SPECIFICATIONS**

Readability = 1 mg

Repeatability = 1 mg

Linearity = ± 0.002 g

Stabilization Time (typical) = 3 s

Tare Range = -520g

Accuracy Class = II

NOTE: exact model may vary, contact sales to confirm current model specifications

**SYRINGE PUMP**  
SYR-PUMP-USB **MOUSE+RAT**

Used for filling bladders at defined rates during cystometry.

The RAZEL™ R100-EC combines the features of our R100-E with the added precision of a computer controlled USB interface. The pump is controlled using MED-CMG software by entering syringe size and desired infusion rate.

- Precision rate control, programmable injections, start/stop & ramp infusion
- LCD display
  - NOTE: To determine flow rates with other syringes, a calibration sheet is included
- Syringe manufacturer & capacity selection
- External & USB operation options
- Accuracy within 2%
- Stop switch disconnects power on syringe depletion
- Infusion rate: 0.6 µL/hr – 467 mL/hr  
0.01 mL/min – 7.78 mL/min





COM-106R2-PC

SG-501

CSG-6080

CSG-6080 cabinet shown with  
IO-115, CANL-430A, IO-206, and IO-200-4

### CYSTOMETRY SYSTEM INTERFACE STARTER

#### CAT-SYST-CYT MOUSE+RAT

Our custom-built data acquisition computer with MED-CMG™ and Cystometry Data Analysis® software installed, comes complete with a 19" wide-screen monitor, speakers, wired keyboard, mouse & pad.

One data acquisition computer can control up to four cystometry systems.

#### SUMMARY OF CONTENTS

- > Pressure transducer calibration kit
- > Computer (w/PCIe card)
- > MED-CMG and Data Analysis Software
- > Tabletop interface cabinet (w/interface card)
- > DC power supply
- > 68-pin cable

### CYSTOMETRY SYSTEM INTERFACE BASIC

#### CAT-SYST-CYT-BASIC MOUSE+RAT

The same contents as the starter interface, only without the computer.



One data acquisition computer can control up to four cystometry systems.

#### COMPUTER REQUIREMENTS

- Windows 7 or newer
- One PCIe slot
- Two USB ports per cystometry station

#### PACKAGE CONTENTS

A = CAT-SYST-CYT		B = CAT-SYST-CYT-BASIC		
#	SKU	DESCRIPTION	A	B
1	CAL-PRESSURE-TRANS	Pressure Transducer Calibration Kit	A	B
1	COM-106-NV	Computer	A	
1	COM-4-PORT-USB-PCI	4 Port USB 2.0 PCIe Card	A	B
1	COM-MONITOR-LCD	Widescreen LED Monitor	A	
1	CSG-6080	Tabletop Interface Cabinet	A	B
1	DIG-744e	PCIe Data Acquisition Card	A	B
1	IO-115	Physiological Interface Card	A	B
1	SG-244-68	Shielded Cable, 68-pin, M/M	A	B
1	SG-501	DC Power Supply 28V, 1A	A	B
1	SOF-551	MED-CMG Cystometry Software	A	B
1	SOF-552	Cystometry Data Analysis Software	A	B



## PRESSURE TRANSDUCER

**TDR-100-1** **MOUSE+RAT**

- Records intra-bladder pressures
- Disposable/reusable



## PRESSURE TRANSDUCER CALIBRATION KIT

**CAL-PRESSURE-TRANS** **MOUSE+RAT**

Inflate the manometer by squeezing the bulb to both reach and hold the desired pressure.

### SUMMARY OF CONTENTS

- > Manometer (reads in units of mmHg)
- > Stopcock
- > Luer connector to easily hook up to the Pressure Transducer (TDR-100-1) for calibration

*Includes conversion table with units of cmH<sub>2</sub>O for reference.*

## INTERFACE CONNECTION CARD

**IO-115** **MOUSE+RAT**


For connecting the interface cabinet to the data acquisition computer via a 68-pin cable.  
Used for recording bladder pressure.



## STANDALONE 28 V DC POWER SUPPLY

**SG-501** **MOUSE+RAT** **120/230 V** **1 A DC**

- Includes two individually switched outputs and two unswitched outputs:
- Switched outputs control accessories that need to be independently turned on and off
- Unswitched outputs power accessories that are always on, like cubicle lights and fans
- Eight (8) individual four-pin Molex® connectors



## COMPUTER

**COM-106R2-PC** **MOUSE+RAT**

## LED MONITOR


**COM-MONITOR-LCD** **MOUSE+RAT**

- Interfaces between software and hardware components to facilitate experimental control and data collection
- Designed to support up to four cystometry packages
- Interface cabinet does not have a built-in 28V DC power supply, sold separately

### SPECIFICATIONS

- Windows 11 (64-bit)
- USB 3.0 (x4), USB 2.0 (x2)
- PCI (x1) / PCIe X16 (x1) / PCIe X1 (x2)
- HDMI (x1) / DisplayPort (x2) / VGA (x1)
- Widescreen monitor w/speakers
- Wired keyboard and mouse





SOF-551

## SMALL ANIMAL CYSTOMETRY “MED-CMG” SOFTWARE


SOF-551 **MOUSE+RAT**

Our powerful data acquisition software for recording cystometrograms (CMGs) from small laboratory animals, such as rats or mice.

- Acquire data from up to four cystometry stations per computer
- Enables control of:
  - bladder filling infusion pump
  - recording bladder pressure
  - pressure transducer calibration
  - annotating event comments during an experiment

### COMPUTER REQUIREMENTS

- Windows 7 or newer
- One PCIe slot
- Two USB ports per cystometry station



SOF-552

## CYSTOMETRY DATA ANALYSIS SOFTWARE

SOF-552 **MOUSE+RAT**

Automatically create a graph of the cystometrogram of the exported data from MED-CMG. Analysis of the data is made possible with readouts of bladder pressure, infusion pump, and scale traces to identify void and non-void contractions.

- X-axis = volume of liquid
- Y-axis = intraluminal pressure of the bladder
- Spikes correspond to the bladder contractions associated with the micturition reflex
- Curve formed by the bottom of the plot reflects the level of pressure necessary to void
- Displays up to four data traces in the main window in an easy-to-view “strip chart”
- Spreadsheets can either be printed or saved
- Calculate key urodynamic parameters:
  - Bladder Pressures
    - Minimum
    - Average
    - Threshold
    - Maximum
  - Inter-micturition Interval
  - Infused Volume

- Voided Volume
- Void Rate Index
- Non-void Contractions
- Isolate & analyze each micturition cycle within a region of interest independently by placing the:
  - Green cursor at the beginning of a cycle,
  - Blue cursor at the threshold of micturition, where the bladder pressure begins to rise sharply
  - Red cursor at the end of the cycle
- Once cursors are in position, the data can exported to a spreadsheet in one click

Data saved to the tab delimited text file includes:

- the voided urine volume
- bladder capacity
- bladder threshold pressure
- minimum bladder pressure during filling
- average bladder pressure prior to micturition
- maximum micturition pressure

### COMPUTER REQUIREMENTS

- Windows 7 or newer





p: 802.527.9724 | f: 802.524.2110 | [www.med-associates.com](http://www.med-associates.com) | e: [sales@med-associates.com](mailto:sales@med-associates.com)