


## WIRELESS RUNNING WHEELS

Low Profile . . . . .	2-3
Vertical. . . . .	4-5
Software . . . . .	6-7





## MONITOR WHEEL RUNNING ACTIVITY WITHOUT RUNNING WIRES.

The current trend of housing mouse colonies in individually ventilated cages (IVCs) has allowed for dramatic increases in housing density, improved environmental control, and better biosecurity. However, the low profile configuration of many mouse IVCs has made it more difficult to add enrichment devices, such as a running wheel, to the mouse home cage. Our wireless low-profile mouse running wheel (ENV-047) represents a state-of-the-art device that enables quantitative measurements of mouse running levels in most standard IVCs. We use the popular Fast Trac™ running wheel surface from BioServ® along with our proprietary wireless transmitter to get running data from the home cage to a computer without any wires or cables.

NOTE: New wheels (ENV-047 & -047V) are compatible with new hub (DIG-807) only. The ENV-047 series wheels will not work with legacy hub (DIG-804).


## LOW-PROFILE WIRELESS RUNNING WHEEL

ENV-047 **MOUSE**

## LOW-PROFILE WIRELESS RUNNING WHEEL **NON-COUNTING**

ENV-044-02 **MOUSE**

- Measure circadian rhythms while adding environmental enrichment
- Quantify running levels to determine how they correlate with behavioral changes
- Fits easily into most standard mouse ventilated rack cages
- Open running surface enables wheel running studies on mice that are instrumented or tethered
- Each hub monitors up to 40 wheels, with four hubs per computer
- Energy efficient and nearly maintenance-free (powered by three AAA batteries)
- Affordable and economical
- Easy to clean: hand wash wheel surface and plastic support base, wipe wheel transmitter with sanitation cloths
- Having spare wheels (ENV-044-01) on hand increases efficiency of cage changes



## MONITOR WHEEL RUNNING ACTIVITY WITHOUT RUNNING WIRES.

Our easy to use & lightweight anodized aluminum wheels fit on any standard wire-topped home cage, and require no tricky set-up or configuration. Perfect for long-term circadian rhythm studies, these wheels are an affordable option for your lab. These wheels are also compatible with our low profile wireless running wheel interface system.

NOTE: New wheels (ENV-047 / -047V) are compatible with new hub (DIG-807) only. The ENV-047 series wheels will not work with legacy hub (DIG-804).


## VERTICAL WIRELESS RUNNING WHEEL

### ENV-047V MOUSE

- Fits onto conventional wire-topped cages
- Measure circadian rhythms with 30 second resolution while adding environmental enrichment
- Enables quantification of running levels and how they correlate with behavioral changes
- Lightweight but heavy duty aluminum vertical wheel revolves on a smooth frame which protrudes through the wire cage top
- Easy-to-clean design
- Designed to be affordable, economical, energy efficient, and nearly maintenance free
- Having spare wheels (ENV-044V-01) on hand increases efficiency of cage changes



ENV-044V-01



## WIRELESS DEVICE USB HUB

### DIG-807 MOUSE

This device connects to the data acquisition computer via a USB cable. The hub must be set to one of eight available addresses and all sensors must be set to the same address as the hub with which they communicate.

- Manage up to 40 wheels per hub, with four hubs per computer
- Hub receives messages once every 30 seconds from the wireless wheels
- Hub relays data to Wheel Manager Software (SOF-860) for data storage and retrieval
- Can be used with both vertical and low-profile wireless running wheels



**Wheel Manager**

File Tools Help

Wheel Sensors

Hub	Id	Name	Time	Count	Battery (v)
1	1	Animal 1	00:00:18:26	37	4.615
1	2	Animal 2	00:00:18:26	12	4.068
1	3	Animal 3	00:00:18:26	24	4.190

Environmental Sensors

Hub	Id	Type	Name	Light (au)	Temperature (C)	Humidity (%)	Time	Battery (v)
1	4	Environmental	Sensor 1	75.000	24.300	39.000	00:00:18:27	4.594

Hubs: 1    Wheels: 3    00:00:18:48    Acquiring

**Export Data**

Source Data File: C:\Documents and Settings\All Users\Documents\20080207\_153527.wls

Start Date: 2/ 7/2008    Time: 3:35:27 PM

End Date: 2/ 7/2008    Time: 3:40:21 PM

Bin Size (min): 1

Output Options

☐ Raw Data only    ☐ Use Sensor Names

Report Format

☒ Date/Time Column    ☒ Include Heading

☐ Time Column    ☒ Include Row Labels

☐ Bin Number Column    ☒ Include Column Labels

☐ Output Nulls as Zero

Wheel Sensors

Hub	Id	Type	Name
<input checked="" type="checkbox"/>	1	0	
<input checked="" type="checkbox"/>	1	1	
<input checked="" type="checkbox"/>	1	2	


Environmental Sensors

Hub	Id	Type	Name
<input type="checkbox"/>	1	4	Humidity (%)
<input type="checkbox"/>	1	4	Temp (C)
<input type="checkbox"/>	1	4	Light (au)

Output File: C:\Documents and Settings\All Users\Documents\20080207\_153527.xls

OK    Cancel

Export data screen with source data selected



Actogram

Periodogram

## WHEEL MANAGER DATA ACQUISITION SOFTWARE

SOF-860 **MOUSE**

- Record running wheel data from wireless running wheels to your computer
- Seamlessly export wheel manager data to Microsoft® Excel or other data analysis software to generate actograms and periodograms
- Export either the entire file or a user-definable time bin/selected time range

## WHEEL ANALYSIS SOFTWARE

SOF-861 **MOUSE**

- Complements our mouse wireless running wheels and Wheel Manager software to analyze circadian rhythms
- Designed for circadian rhythm analysis in real-time during acquisition or offline from saved Wheel Manager data
- Generates actograms and periodograms



p: 802.527.9724 | f: 802.524.2110 | [www.med-associates.com](http://www.med-associates.com) | e: [sales@med-associates.com](mailto:sales@med-associates.com)