




VeriFEED Dispenser

ENV-204

USER MANUAL
DOC- 337
Rev.1.0



Copyright ©2023 All Rights Reserved

Phone: 802.527.2343

Fax: 802.524.2110

Med Associates Inc.
P.O. Box 319
St. Albans, VT 05478
www.med-associates.com

CONTENTS

CHAPTER 1 INTRODUCTION	1
Removable Hopper	1
Feeder Bracket.....	1
Exchangeable Pellet Discs and Improved Pellet Handling	1
Changing of the Guard	1
Figure 1.1 Connections and Controls	2
Choice of Mounting Options	3
Figure 1.2 Mounting and Accessories.....	3
CHAPTER 2 GETTING SET-UP	4
Assembling VeriFEED connections.....	4
Figure 2.1 Assembling VeriFEED	4
MedINPUT	4
Figure 2.2 The pinout of the 3-pin MedOPERATE and MedINPUT Molex connector	4
Figure 2.3 Hopper Removal Cleaning Assembly and Attachment.....	5
Minimum Fill Line.....	7
Figure 2.4 Image of Minimum Fill line	7
Receptacle Hose preparation.....	7
Figure 2.5 Image of Receptacle Hose preparation	7
MedOPERATE.....	7
Feeder Empty Alarm.....	8
Confirming operation in MedTest and in place	8
Hopper Assembly and Attachment	8
Getting Started	8

CHAPTER 3 SAMPLE MED-PC MSN PROGRAMS	9
K-pulse (#k).....	9
Figure 3.1 Send K, R or Start Signal Button	9
Figure 3.2 Send Signals to Boxes Window	9
Examples and Descriptions	10
Example 3.1 Single Pellet Delivery	11
Example 3.2 Multiple Pellet Delivery	12
Example 3.3 Single Pellet with Counter	13
Example 3.4 Multiple Pellet with Counter	14
Example 3.5 Multiple Pellet with Counter and Verification.....	15
CHAPTER 4 SAFETY AND FAIL-SAFES	16
Cleaning	16
Warranty & Repair	16
Mechanical Maintenance.....	16
Table 4.1 Mechanical Properties Table	17
Contact Information.....	17

CHAPTER 1 INTRODUCTION

Med Associates recommends reading this manual prior to operating the VeriFEED Dispenser.

The VeriFEED is a landmark design in pellet dispensers. Reward variability in experiments should be based on experimental design not on performance of equipment. Design features introduced in this model ensure reliable verified repeatable delivery and easy maintenance.

Removable Hopper

The pellet hopper snaps to the *feeder bracket* enclosure with strong magnets and opens into two *hopper halves* with the flip of a *clip*. The upright angled fill mouth and simple *lid* make it easy to see, empty, clean, fill, and weigh.

Feeder Bracket

The Feeder Bracket *stepper motor* has improved torque control, and an Infrared sentry with smart gain technology detects pellets even as dust accumulates. *IR Sentry* and stepper motor protocols coordinate action to ensure delivery, and self-clear most blockages. The VeriFEED sends a verification signal to *Med-PC* after each pellet is counted. In the rare event that it is unable to dispense a pellet, it will signal with a flashing light the user can also have *Med-PC* notify when no pellet is dropped.

Exchangeable Pellet Discs and Improved Pellet Handling

The precision tooled *ENV-204D-20* (20 mg) or *ENV-204D-45* (45 mg) pellet discs for grain or sucrose provide a gentle profile for handling pellets. The angle of operation reduces interference which could cause grinding together or crumbling of pellets. Pellet discs are sold separately.

Changing of the Guard

The VeriFEED replaces all of our previous 20mg and 45mg dispenser lines.

This includes:

ENV-203-20	Pedestal Dispenser for 20mg Pellets
ENV-203-20IR	Pedestal Dispenser with IR Sentry for 20mg Pellets
ENV-203-45	Pedestal Dispenser for 45mg Pellets
ENV-203-45IR	Pedestal Dispenser with IR Sentry for 45mg Pellets
ENV-203M-20	Modular Dispenser for 20 mg Pellets
ENV-203M-45	Modular Dispenser for 45mg Pellets
ENV-203M-45IR	Modular Dispenser with IR Sentry for 45mg Pellets

Figure 1.1 Connections and Controls

The VeriFEED Dispenser

- A = Hopper
- B = Status Light
- C = Bracket
- D = Power and Control Port
- E = Post
- F = Hopper Lid
- G = Pellet Drop Tube
- H = Operate Button




Choice of Mounting Options

The VeriFEED comes ready to be mounted on:

- Any Med Associates modular panel including Classic Mouse, Wide Mouse, or Standard Rat.
 - A 90° mount that attaches to the same modular panel but turns the VeriFEED hopper to face the front of the sound attenuating cubicle (SAC).
 - A pedestal stand that sets it off to the side and can attach to the environment base.
- NOTE: Mounting options are sold separately.

Figure 1.2 Mounting and Accessories



ENV-204P

Pedestal Mount



ENV-204M90M

Modular Chamber Panel Mount w/ 90° turn for Classic Mouse

ENV-204M90WM

Modular Chamber Panel Mount w/ 90° turn for Wide Mouse

ENV-204M90R

Modular Chamber Panel Mount w/ 90° turn for Standard Rat



ENV-204MM

Modular Chamber Panel Mount for Classic Mouse

ENV-204MWM

Modular Chamber Panel Mount for Wide Mouse




ENV-204MR

Modular Chamber Panel Mount for Standard Rat

CHAPTER 2 GETTING SET-UP

Assembling VeriFEED connections

Figure 2.1 Assembling VeriFEED



The VeriFEED is designed to be controlled through a MedConnect style connection panel such as the OC-112 or the SG-716B. See the *Omni Control Interface Manual* or the *Smart Control Interface Manual* for more information on that operation. In brief, the power and control port can be found at the lowest part of the bracket ("D" in Figure 1.1 Connections and controls). The SG-224A cable separates into two connections, a MedINPUT and a MedOPERATE. Both component connections to the VeriFEED terminate in a 3 pin Molex style connector.


MedINPUT



The MedINPUT connector carries simple signals to Med-PC. The VeriFEED dispenser sends a 20 ms pulse to the connection panel when a pellet passes the IR sentry through the pellet drop tube. This can then be used to trigger reaction events in *Med State Notation protocols (MSN)* or notify the user that a feeder is empty (See CHAPTER 3 SAMPLE MED-PC MSN PROGRAMS). The feeder will operate normally if the MedINPUT line is not connected but will not send any responses to Med-PC.

Figure 2.2 The pinout of the 3-pin MedOPERATE and MedINPUT Molex connector

	Control Panel connections	VeriFEED connection
MedOPERATE		
1 Output Ground (28 Volt to Ground)	Green	1
2 Operate Active low (+28 Volts)	Blue	2
3 Unit Power +28 Volts DC	Orange	3
MedINPUT (VeriFEED Pellet signal)		
4 Input Ground	Black	4
5 Signal Line Active low (+28 Volts)	White	5
6 +28 Volts DC (Connection only)	Red	6

Figure 2.3 Hopper Removal Cleaning Assembly and Attachment





Minimum Fill Line

When the pellet dispenser runs close to empty, some pellet geometry can work against the operation of the machine to create matrixed or linked pellets that lift away from the entry points. The VeriFEED will still dispense pellets with the same accuracy as before, but the time between dispenses may slow as the feeder works to counter this interference. The feeder will run most efficiently when there is some weight of pellets above the entry to the disc portion of the hopper. The minimum fill line on the hopper shows where this point begins, and where less efficient dispensing may start to happen.


Figure 2.4 Image of Minimum Fill line



Receptacle Hose preparation

A flexible hose line provides a seamless connection between the drop point and the pellet receptacle. It is cut to size to give a secure and un-kinked path to the pellet receptacle. If it becomes necessary to replace this hose the new hose can grip very firmly to the pellet drop tube. We suggest the operator cut a single slit up the last 2-3 cm of line the pellet drop tube will insert into. To remove the hose a firm grip at the base of drop tube and a rocking twist should remove the hose fairly easily. Note when attaching the hose that the pellet drop tube can be lifted upward. Press it down so that the retaining clip comes into contact with the housing and therefore is seated completely into position.

Figure 2.5 Image of Receptacle Hose preparation



MedOPERATE

The MedOPERATE is simultaneously the power for the VeriFEED and the signal for the device to operate.

Med-PC sends active-low operate signals. Meaning that the signal is active when the voltage drops to ground on the operate line. A pulse to operate is produced between the operate signal starting and stopping, meaning power returns to the operate pin from ground. This reduces latency in operation, and helps confirm connection integrity, but care must be taken if damage is done to the operate pin (2) as a short between this line and the Unit Power pin (3) will trigger the feeder.

The VeriFEED is capable of retaining a very large number of sent signals to operate in queue while the feeder dispenses. In laboratory tests it parses signals up to 20 times a second. Barring interference, the feeder will dispense a pellet within one second the majority of the time.

Feeder Empty Alarm

If the pellet dispenser has completed its logic sequence to deliver a pellet, and is still unsuccessful, the status light flashes (See Figure 1.1 Connections and controls p.2). See SAMPLE MED-PC MSN PROGRAMS for information on using this feature to warn you of an empty hopper using your Med-PC protocol.

Confirming operation in MedTest and in place

When the VeriFEED feeder has been connected to the ports the user intends, then the operation can be checked as soon as the chamber is powered by pressing the operation button on the bracket bottom (Figure 1.1 Connections and controls). If the hopper is filled, it will deliver a pellet with each press of the button. You can also stimulate the feeder through MedTest by opening the control panel for the chamber it is attached to and double clicking the port number

Hopper Assembly and Attachment

The Hopper is designed to be intuitively assembled without tools. Insert the curved tabs into the hinge posts and line up the ridges and grooves along the perimeter of the hopper. Insert the pellet disc into the disc housing on either side of the hopper with the spacer post aiming up toward the pellet reservoir and the pellet entry holes against the bottom of the hopper. Close the two halves together over the pellet disc. Bring the clip over the screw post and press the tension lever of the pressure clip against the housing to snap into place. To attach the hopper, bring the hopper toward the feeder. Visually align the pellet entry point with the pellet drop tube and the square disc hub bore with the square axle drive of the motor. It may take some movement to find the setting, as the disc rotates freely within the hopper till seated. Slide in till the magnets catch and hold the hopper in place. You can now fill the hopper with appropriately sized pellets and cover the top with the included cap.

Getting Started

Once the feeder is assembled, connected, and filled as described above, it is ready to be used as a verified reliable pellet dispenser in your behavior protocol.


CHAPTER 3 SAMPLE MED-PC MSN PROGRAMS

Refer to the Med-PC User's Manual and Med-PC Programmer's Manual for more detailed information and examples for coding Med State Notation protocols.

K-pulse (#k)


The following five MedState Notation example programs use a "K1 pulse" issued via the "Send Signals to Boxes" dialog. To launch the Send Signals to Boxes dialog, click the "Traffic Light" toolbar icon, or open the "Sessions" menu and select "Signals (START, #K, #R)...".

Figure 3.1 Send K, R or Start Signal Button



Choose the "K pulse" Signal option, select the Box running the example, and click "Issue" to deliver a K1 pulse.

Figure 3.2 Send Signals to Boxes Window



A K1 pulse is excellent for supplying inputs when testing MedState Notation programs in MED-PC without relying on using actual test chamber devices. When the procedure is used in a real scenario, a lever or other input device will replace the K1 pulse.

In each example the ENV-204's Output Line is connected to Output 1, and the Input Line is connected to Input 1 on a SmartCtrl connection panel. The Output Line triggers the operation of the pellet delivery mechanism, and the Input Line indicates the pellet feeder's Infrared pellet detector sensed an object dropping from the feeder.

Examples and Descriptions

The examples are presented in order of increasing complexity like the Tutorials provided with the MED-PC Trans programming utility.

Example 1 (3.1) delivers a single-pellet reward. There is no counting of pellets delivered and pellet delivery verification (no error-checking).

Example 2 (3.2) delivers a multi-pellet reward. As in the first example, there is no counting of pellets delivered nor delivery verification performed.

Example 3 (3.3) delivers a single-pellet and counts the number of pellets delivered by adding a dedicated pellet-counting State Set. This added State Set uses the input sent from the ENV-204 when its Infrared detector "sees" a pellet falling through the delivery chute.

Example 4 (3.4) delivers a multi-pellet reward and displays the number of rewards and pellets delivered.

Example 5 (3.5) delivers a multi-pellet reward, display the number of rewards and pellets delivered, and adds error-checking and reports undelivered pellets.

Example 3.1 Single Pellet Delivery

```

1  \Example 1: deliver a single-pellet reward with each K1 pulse
2
3  \OUTPUTS
4  ^VeriFEED    = 1
5  \ENV-204 VeriFEED connected to Output 1 on DIG-716B Connection Panel
6
7  S.S.1,
8  \State Set 1
9  S1,
10     0.01": SHOW 1, ISSUE K1 to deliver pellet, 1 ---> S2
11     \Instructions for use: issue K1
12 S2,
13     #K1: ON ^VeriFEED ---> S3
14     \When K1 pulse received, turn on VeriFEED output
15 S3,
16     0.02": OFF ^VeriFEED ---> S2
17     \Turn off VeriFEED output and go back to S2 to look for next K1 pulse

```

Example 3.1 shows the simplest method of using the ENV-204 pellet feeder: turn on the output connected to the feeder, and then turn that output off, repeat as desired. When using this program, experiment by issuing K1 pulses in rapid succession by clicking the “Issue” button on the “Send Signals to Boxes” dialog. You will be able to issue commands to the feeder faster than it can mechanically deliver the pellets. Notice the feeder will keep track of the number of commands it received, and how many pellets it delivered, and perform its delivery protocol until those pellets are delivered. If the pellets are not delivered in about 80 seconds, the feeder enters the “empty hopper” state and flashes its error-indicating LED. This error state is not handled in this example but is managed in Example 3.5.

Example 3.2 Multiple Pellet Delivery

```

1  \Example 2: deliver a multi-pellet reward with each K1 pulse
2
3  \OUTPUTS
4  ^VeriFEED = 1
5  \ENV-204 VeriFEED connected to Ouput 1 on DIG-716B Connection Panel
6
7  \CONSTANTS
8  ^RewardSize = 5  \Deliver 5 pellets per reward
9
10 S.S.1,
11 \State Set 1
12 S1,
13   0.01": SHOW 1, ISSUE K1 to deliver pellets, ^RewardSize ----> S2
14   \Instructions for use: issue K1
15 S2,
16   #K1: SET P = 0 ----> S3
17   \When K1 pulse received, set counter P to 0
18 S3,
19   0.02": IF P < ^RewardSize [@t, @f]
20   \If P less than RewardSize (5 in this example)
21   @t: ADD P; ON ^VeriFEED ----> S4 \Then increase P and deliver pellet
22   @f: ----> S2 \Else wait for next K1 pulse
23 S4,
24   0.01": OFF ^VeriFEED ----> S3 \Turn off VeriFEED output

```

Example 3.2 builds upon the previous example by issuing multiple pellets (five) with each K1 pulse. An "If Statement" in State 3 controls the logic and counting of pellets delivered. The code does not look for inputs from the ENV-204 indicating a pellet was physically delivered, but rather uses an internal counter "P" to keep track of pellets delivered.

Example 3.3 Single Pellet with Counter

```

1  \Example 3: deliver a single-pellet reward with each
2  \           K1 pulse and count the number of pellets delivered
3
4  \OUTPUTS
5  ^VeriFEED = 2
6  \ENV-204 VeriFEED connected to Output 1 on DIG-716B Connection Panel
7
8  \INPUTS
9  ^Delivered = 1
10 \ENV-204 VeriFEED read connected to Input 1 on DIG-716B Connection Panel
11
12 \VARIABLES
13 \D - Delivered pellets detected by feeder IR detector
14 S.S.1,
15 S1,
16   0.01": SHOW 1, ISSUE K1 to deliver pellet, 1 ----> S2
17   \Instructions for use: issue K1
18 S2,
19   #K1: ON ^VeriFEED ----> S3
20   \Turn ON VeriFEED output to dispense pellet
21 S3,
22   0.02": OFF ^VeriFEED ----> S2
23   \Turn OFF VeriFEED output to be ready for next delivery
24
25 S.S.2,
26 \State Set 2 looks for an input from the feeder's IR detection line
27 S1,
28   #R^Delivered: ADD D; SHOW 2, Delivered Pellets, D ----> S1
29   \Add to counter, show pellets detected

```

Example 3.3 uses the same State Set 1 as Example 3.1 to deliver a single pellet with each K1 pulse. A constant for the ENV-204 input connector is added on line 9 above. Also, State Set 2 (lines 25-29) is added to monitor the ENV-204 Med Input line. When an input is detected in State Set 2, State 1, the VeriFEED has “seen” a pellet dropping from the dispenser into the pellet delivery tube. SHOW field 1 indicates the number of pellets seen leaving the ENV-204 feeder instead of relying on an internal counter as in Example 3.2.

Example 3.4 Multiple Pellet with Counter

```

1  \Example 4: deliver a multi-pellet reward with each
2  \           K1 pulse and count number of pellets delivered
3
4  \OUTPUTS
5  ^VeriFEED = 1
6  \ENV-204 VeriFEED connected to Output 1 on DIG-716B Connection Panel
7  \INPUTS
8  ^Delivered = 1
9  \ENV-204 VeriFEED read connected to Input 1 on DIG-716B Connection Panel
10 \VARIABLES
11 \D - Delivered pellets detected by feeder IR detector
12 \P - Pellets per Reward
13 \R - Rewards issued (5 pellets per Reward)
14 \CONSTANTS
15 ^RewardSize = 5    \ deliver 5 pellets per reward
16
17 S.S.1,
18 S1,
19     0.01": SHOW 1, ISSUE K1 to Reward pellets, ^RewardSize ---> S2
20     \ Instructions for use: issue K1
21 S2,
22     #K1: SET P = 0 ---> S3
23 S3,
24     0.02": IF P < ^RewardSize [@t, @f]
25     \ If counter is less than pellets to be delivered...
26         @t: ADD P; ON ^VeriFEED ---> S4
27         \Then increase counter and turn on VeriFEED output
28         @f: ---> S2
29         \Else wait for next K1 pulse
30 S4,
31     0.01": OFF ^VeriFEED ---> S3
32     \Turn off VeriFEED output to be ready for next delivery
33
34 S.S.2,
35 \State Set 2 looks for an input from the feeder's IR detection line
36 S1,
37     #R^Delivered: ADD D; SHOW 2, Pellets Delivered, D ---> S1
38     \Show number of pellets detected by IR detector
39
40 S.S.3,
41 \State Set 3 counts the times a multi-pellet reward was requested
42 S1,
43     #K1: ADD R; SHOW 3, Rewards Requested, R ---> S1
44     \Show number of rewards requested

```

Example 3.4 builds upon Example 3.2 by adding a pellet-delivered counting routine in State Set 2 and reward-requested counter in State Set 3. Five pellets are delivered with each K1 pulse, the number of K1 pulses is displayed in SHOW field 3, and the number of pellets dispensed is printed in SHOW field 2.

Example 3.5 Multiple Pellet with Counter and Verification

```

1  \Example 5: delivers a multi-pellet reward with each K1 pulse, counts the delivered
2  \      pellets, and reports undelivered pellets
3
4  \OUTPUTS
5  ^VeriFEED = 1
6  \ENV-204 VeriFEED connected to Output 1 on DIG-716B Connection Panel
7
8  \INPUTS
9  ^Delivered = 1
10 \ENV-204 VeriFEED read connected to Input 1 on DIG-716B Connection Panel
11
12 \VARIABLES
13 \D - Delivered pellets detected by feeder IR detector
14 \P - Counter of Pellets per Reward
15 \R - Rewards issued (5 pellets per Reward)
16
17 \CONSTANTS
18 ^RewardSize = 5  \Number of pellets to deliver (5 in this example)
19
20 S.S.1,
21 \State Set 1 controls pellet delivery logic
22 S1,
23   0.01": SHOW 1, ISSUE K1 to Reward pellets, ^RewardSize ---> S2
24   \Instructions for use: issue K1
25 S2,
26   #K1: SET P = 0; ADD R; SHOW 3, Rewards Requested, R ---> S3
27   \Look for the K1 pulse, reset P counter
28 S3,
29   0.01": IF P < ^RewardSize [@t, @f]
30   \If counter is less than pellets to be delivered..
31   @t: ADD P; ON ^VeriFEED -> S4
32   \Then increase counter and turn on VeriFEED output
33   @f: ---> S2  \else wait for the next K1 pulse
34 S4,
35   0.02": OFF ^VeriFEED -> S5
36   \Turn off VeriFEED output to be ready for next delivery
37 S5,
38   #R^Delivered: ADD D; SHOW 2, Pellets Delivered, D -> S3
39   \If input detected from IR detector, increase counter
40   80": ADD X; SHOW 4, ERROR, X -> S2
41   \No input after 80 seconds: increase error counter, wait for K1

```

Example 3.5 puts all the previous examples together and delivers a multi-pellet reward with each K1 pulse, SHOWs the number of pellets physically delivered, and displays a warning of pellet non-delivery. The ENV-204 will perform a sequence of movements (forwards, backwards, and vibration) to deliver a pellet. The ENV-204 pellet delivery sequence lasts approximately 80 seconds. Example 5, State 5 looks for a pellet delivery for up to 80 seconds. If no pellet is seen after this duration, an error counter "X" is incremented and displayed in SHOW field 4. SHOW field 2 displays the pellets delivered, SHOW field 3 displays the number of rewards requested.

CHAPTER 4 SAFETY AND FAIL-SAFES

Cleaning

Do not submerge any portion of the feeder in liquid. Clean feeder parts using a soft cloth that is dampened with water and/or detergent. Clean the pellet drop tube by passing a dry cotton swab through it to remove pellet dust. The feeder bracket enclosure is stainless steel and aluminum. Modular Panels and mounts are stainless steel. The pellet disc and hopper are injection molded polycarbonate with strong nickel-plated neodymium magnets and brass heat insets for firm attachment of screws. Magnets are vulnerable to degradation due to heat and moisture. Should fluids enter the feeder internals, disconnect the power, and invert the feeder to drain. Set the feeder in a warm dry area to air dry thoroughly before restoring power. It is recommended to empty and wipe out the hopper at the end of each use period (See Hopper Assembly and Attachment p.8) to improve the quality of the pellets and how well they are dispensed.



DO NOT AUTOCLAVE



DO NOT SUBMERGE

Warranty & Repair

This device is protected under a limited manufacturer warranty for two years from the date of purchase.

Med Associates continues to support and repair products outside the warranty period indefinitely, for a fee, if we are logistically able. We are proud to continue to support instruments manufactured in our opening production line from 1972.

If you are interested in further information about this device, its component parts, accessories, or integrating it within your application, please contact Med Associates Support (See p.17).

Mechanical Maintenance

High quality pellets and clean pellet preparations improve feeder performance. It is recommended that the hopper be removed, emptied, opened, and wiped out with a clean dry cloth between uses. Parts should be inspected for wear; a worn part should be replaced by Med Associates (Contact Med Associates support see p.17 below).

Table 4.1 Mechanical Properties Table

	WEIGHT*	DIMENSIONS l* x w* x h*	MATERIALS (primarily)
Hopper Left	3 oz / 84 g	4.4 in x 2.2 in x 4.4 in 11 cm x 5.5 cm x 11 cm	Polycarbonate
Hopper Right	3.2 oz / 90 g	5.6 in x 2.2 in x 4.4 in 14(11) cm x 5.5 cm x 11 cm (w/o clip)	Polycarbonate, Aluminum
Hopper Lid	1.7 oz / 48 g	4.2 in x 3 in x 1 in 10.5 cm x 7.5 cm x 2.5 cm	Polycarbonate
Pellet Disc	1.4 oz / 40 g	3.6 in x 3.6 in x 1 in 9 cm x 9 cm x 2.5 cm	Polycarbonate
Total Hopper	8.5 oz / 242 g	4.4 in x 4.4 in x 4.4 in 11 cm x 11 cm x 11 cm	Polycarbonate
Bracket	17.2 oz / 488 g	6.7 in x 5.1(1.6) in x 2.9 in 17 cm x 13(4) cm x 7.5 cm (w/o post)	Steel, Aluminum
VeriFEED (no mount)	25.7oz / 730 g	7.7 in x 4.4 in x 7.5 in 19.5 cm x 11 cm x 19 cm	Steel, Aluminum, Polycarbonate

*Measures are approximated for convenience in space planning and subject to change

MAX DELIVERY RATE 20 mg /min	150 pellets
Max Delivery Rate 45 mg /min	120 pellets
Hopper Capacity 20 mg	4500 pellets
Hopper Capacity 45 mg	3000 pellets
VeriFEED Delivery Accuracy Rate	>99.9%**

**In continuing trials of 1000 expected deliveries, less than 1 delivery is missed.

Contact Information

Please contact Med Associates, Inc. for information regarding any of our products. Med Associates Support is available to answer technical questions weekdays between 8:00am and 4:30pm ET by phone at 802-527-2343, or email at support@med-associates.com. We look forward to collaborating with you.

Visit our website at www.med-associates.com
For pricing inquiries, email sales@med-associates.com

INDEX

2
20 ms pulse4

3
3 pin Molex style connector ...4

A
Active low signal7

B
Bracket2

C
clip1
connection panel4

E
ENV-203-201
ENV-203-20IR1
ENV-203-451
ENV-203-45IR1
ENV-203M-201
ENV-203M-451
ENV-203M-45IR1
ENV-204D-201

NOTES:

ENV-204D-451
Environment base3
Exchangeable Pellet Disc1

F
Feeder Bracket1

H
Hopper1, 2
hopper halves1
Hopper Lid2

I
infrared sentry1
IR sentry4
IR Sentry1

L
lid1

M
Med State Notation (MSN)4
MedConnect4
MedINPUT4
Med-PC1, 4
Minimum Signal Resolution7

Modular Panel3
Mounting Options3

O
OC-1124
Operate Button2
OPERATE pin7
OPERATE signal7

P
Pedestal Mount3
Pellet Drop Tube2
Post2
Power and Control Port2

S
SG-224A cable4, 7
SG-716B4
Sound attenuating cubicle (SAC)3
Status Light2
stepper motor1

V
VeriFEED1
Voltage pin7

Copyright ©2023 All Rights
Reserved

Phone: 802.527.2343

Fax: 802.524.2110

Med Associates Inc.
P.O. Box 319
St. Albans, VT 05478
www.med-associates.com

This Manual was prepared by:
Med Associates Staff
With special thanks to the contributions of

Seth Truran
Application Scientist

Kevin Deutermann
Graphic Designer

Connor Harrigan
Mechanical Engineer

Richard Wheeler
Electrical Engineer


Gary Bamberger
Software Engineer

Andy Lebeis
Software Engineer

Jon Beyor
Support Manager

Bobbi Bushey
Sales Manager

Chris Johnson
Engineering Manager



Copyright ©2023 All Rights Reserved

Phone: 802.527.2343

Fax: 802.524.2110

Med Associates Inc.
P.O. Box 319
St. Albans, VT 05478
www.med-associates.com