


ULTRA SENSITIVE MOUSE LEVER

ENV-310M

HARDWARE SHEET

Molex Pin #	Wire Color	Function
1	Black	28 Volt Ground
2	White	Lever Response (Input)
3	Red	+28 Volts



Connect the Response cable to any available Input on a Standard MED Connection Panel.

DOC-154

Rev. 1.2

Copyright © 2012
All Rights Reserved

Med Associates Inc.
P.O. Box 319
St. Albans, Vermont 05478

Phone: 802.527.2343
Fax: 802.527.5095
www.med-associates.com