

# **FIXED WATER MAZE PLATFORM**

ENV-596

OPERATING INSTRUCTIONS

**DOC-191**

**Rev. 1.1**

Copyright ©2014

All Rights Reserved

Med Associates Inc.

P.O. Box 319

St. Albans, Vermont 05478

Phone: 802.527.2343

Fax: 802.527.5095

[www.med-associates.com](http://www.med-associates.com)

## CHAPTER 1 | OPERATION

### Raising the Platform

1. Turn the T-screws counterclockwise to loosen.
2. Raise the platform to the desired height.
3. Turn the T-screws clockwise to tighten.



### Lowering the Platform

1. Turn the T-screws counterclockwise to loosen.
2. Lower the platform to the desired height.
3. Turn the T-screws clockwise to tighten.

## APPENDIX A | CONTACT INFORMATION

Please contact MED Associates, Inc. for information regarding any of our products.

Visit our website at [www.med-associates.com](http://www.med-associates.com) for contact information.

For technical questions, email [support@med-associates.com](mailto:support@med-associates.com).