

# **TWEETER ASSEMBLY**

ENV-324DM

ASSEMBLY MANUAL

**DOC-275**

**Rev. 1.0**


Copyright ©2015  
All Rights Reserved

Med Associates Inc.  
P.O. Box 319  
St. Albans, Vermont 05478

Phone: 802.527.2343  
Fax: 802.527.5095  
[www.med-associates.com](http://www.med-associates.com)

## ASSEMBLY INSTRUCTIONS

- Remove the four #6 screws from the top of the ENV-307A mouse chamber.
- Replace the existing top with the ENV-324DM tweeter assembly as shown below.
- Replace the four #6 screws.



## CONTACT INFORMATION

Please contact Med Associates, Inc. for information regarding any of our products.

Visit our website at [www.med-associates.com](http://www.med-associates.com) for contact information.

For technical questions, email [support@med-associates.com](mailto:support@med-associates.com).