

SOUND ATTENUATING CUBICLE (SAC) WITH VIEWING WINDOW

ENV-018MD-W

PRODUCT MANUAL

DOC-280

Rev 1.0

Copyright ©2013
All Rights Reserved

Med Associates Inc.
P.O. Box 319
St. Albans, Vermont 05478

Phone: 802.527.2343
Fax: 802.527.5095
www.med-associates.com

notes


Table of Contents

| | |
|--|----------|
| Chapter 1 Introduction..... | 1 |
| Appendix Contact Information..... | 2 |

CHAPTER 1 | INTRODUCTION

The Sound Attenuating Cubicle (SAC) with a clear polycarbonate viewing window in the left door is ideal for group viewing or videotaping applications. The viewing window allows researchers to observe an animal in the cubicle without opening the door. The SAC is manufactured using Multi Density Fiberboard (MDF), a high density wood composite with thermally infused high pressure laminate inside and out. A factory installed 28V DC fan is included.

Figure 1.1 - Sound Attenuating Cubicle (SAC) with Viewing Window



APPENDIX | CONTACT INFORMATION

Please contact MED Associates, Inc. for information regarding any of our products.

Visit our website at www.med-associates.com for contact information.

For technical questions, email support@med-associates.com.