

CYSTOMETRY FITTINGS KIT FOR RAT

CYT-R-FIT-KIT

This kit contains the following:

Quantity	Catamount Part Number	Description
1 pkg (10)	TUB-NDL-18G-S	18 Gauge Needle Stub
1 pkg (10)	TUB-NDL-23G-S	23 Gauge Needle Stub
1 pkg (10)	TUB-CYT-CATH-R	Catheters (see Disclaimer below)
1 pkg (24")	TUB-TYGON-040	.040" ID Tygon Tubing
1 pkg (10)	LUER-SLIP-COUPLER-MALE-PK	Luer Slip Coupler
1 pkg (10)	LUER-STOPCOCK-3WY-PK	3-Way Stopcock
1 pkg (5)	SYR-20	20 ml (cc) Syringe

DISCLAIMER

CYSTOMETRY RAT CATHETERS

This kit contains 10 rat bladder catheters.

Contents are not sterile. During the manufacture process, catheters are flushed with a small amount of ethyl alcohol. Thus, catheter lumens may contain ethyl alcohol. Catheters should be flushed with bladder filling solution prior to use. In most cases, catheters should be sanitized/sterilized prior to implantation. Consult with your Institution Animal Care and Use Committee prior to implanting these catheters into animals.

DOC-162
Rev. 1.4

Copyright © 2024
All Rights Reserved

Med Associates, Inc.
166 Industrial Park Road
Fairfax, VT 05454
www.med-associates.com