

NEW FROM **med associates inc**
PROJECTED RELEASE Q1 2027 REQUEST A QUOTE TODAY!



- Capture analysis-ready footage from up to 16 cameras with one PC
 - With distributed computing, use any PC running Windows 10 or 11
- Tie video to behavior data with Med-PC, and sync timecode to Med-PC's common clock
- Match data events and trials to specific video segments with synchronized timecodes
 - Two separate counters can be used for time stamping
- File everything under the correct trial, segment, and date with user-specified naming convention
- Annotate the video file using Med-PC
 - Can be overlaid on the video or saved in a separate annotation file
- Start and stop recording via Med-PC or Med Video application
- Med Video hardware can be used outside of a Med-PC environment to work with a wide range of hardware (ex. TTL triggers)

**>CONTACT OUR SALES TEAM TO
EXPLORE HOW MED VIDEO CAN
BE ADDED TO YOUR SETUP!**

	MED VIDEO
COLOR	Monochrome
FORMAT	H.264 / MPEG-4 AVC (*.mp4)
FRAME RATE*	Up to 60 FPS
RESOLUTION*	Up to 1080p
TIMECODE	Millisecond (0.000s)
LIGHTING	White or Near Infrared (NIR)
CONTINGENCY STORAGE	24 hrs @ highest quality on built-in microSD
CAMERA OPTIONS	VID-CAM-DVC (720p) VID-CAM-MONO-5 (720p) VID-CAM-MONO-6 (720p) VID-CAM-MONO-7 (1080p)

*NOTE: Varies by camera and settings